
Laser μ -FabTM

Operation Manual

2024.06.05

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5)	DEVICE CONTROL	45
1)	SCIENTIFIC CHART	46
6)	MEASUREMENT	48
7)	CAD GENERATOR	49
7)	IMAGE TO DXF	52
8)	ALIGNMENT TOOL – MOTION CONTROL	53
8.	참고자료	60
1)	TROUBLE SHOOTING	60
9.	레이저 시스템 도면	66

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▶ IEC (International Electrotechnical Commission)

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▶ ANSI(American National Standards Institute)

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▶ FDA/CDRH(Center for Devices and Radiological Health)

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ANSI and IEC Laser classification

ANSI and IEC laser classification	Class 1		Class 2	
	Class 1	Class 1M	Class 2	Class 2M
sub-class				
U.S. FDA laser classification	Class I	No special FDA class	Class II	No special FDA class
Human-accessible laser power (for visible light)	For visible light, emits beam less than 0.39milliwatts, or beam of any power is inside device and is not accessible during operation.		Emits visible beam of less than 1 milliwatts.	
Caution/warning indication	No special caution/warning indication		No special caution/warning indication	
Label descriptive text		DO NOT VIEW DIRECTLY WITH OPTICAL INSTRUMENTS	DO NOT STARE INTO BEAM	DO NOT STARE INTO BEAM OR EXPOSE USER OF TELESCOPIC OPTICS

ANSI and IEC laser classification	Class 3		Class 4
	Class 3R	Class 3B	Class 4
sub-class			
U.S. FDA laser classification	Class IIIa (definition is different but results are similar)	Class IIIb	Class IV
Human-accessible laser power (for visible light)	For visible light, emits beam between 1 and 4.99 milliwatts	For visible light, emits beam between Class 3R limit(e.g.5milliwatts)and 499.9milliwatts	For visible light, emits beam of 500 milliwatts (1/2 Watts) or more
Caution/warning indication	<CAUTION>	<WARNING>	<DANGER>
Label descriptive text	AVOID DIRECT EYE EXPOSURE	AVOID ECPOSURE TO BEAM	AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION

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


 <p>The label consists of a triangular warning symbol with a black border and a black sunburst icon in the center. Below the triangle is a rectangular box with a black border containing the text: "LASER RADIATION", "DO NOT STARE INTO BEAM", and "CLASS 2 LASER".</p>	<p>Class 2 등급의 레이저이며 광선을 직접적으로 보지 말 것을 경고하는 표시.</p>
 <p>The label is rectangular with a black border. At the top is a red oval with the word "DANGER" in white. Below this is a red sunburst icon. To the right of the icon, the text reads: "LASER RADIATION - AVOID DIRECT EYE EXPOSURE", "(LASER TYPE, WAVELENGTH, PULSE DURATION, MAXIMUM OUTPUT)", and "CLASS 3a LASER PRODUCT".</p>	<p>Class 3a(R) 등급의 레이저이며 직접적으로 눈에 노출되면 위험함을 표시.</p>
 <p>The label is rectangular with a black border. At the top is a red oval with the word "DANGER" in white. Below this is a red sunburst icon. To the right of the icon, the text reads: "LASER RADIATION AVOID DIRECT EXPOSURE TO BEAM" and "CLASS 3b LASER".</p>	<p>Class 3b 등급의 레이저이며 직접적으로 눈에 노출되면 위험함을 표시.</p>
 <p>The label consists of a triangular warning symbol with a black border and a black sunburst icon in the center. Below the triangle is a rectangular box with a black border containing the text: "LASER RADIATION", "AVOID EXPOSURE TO BEAM", and "CLASS 3B LASER".</p>	<p>Class 3b 등급의 레이저이며 직접적으로 눈에 노출되면 위험함을 표시.</p>

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Safety Instructions

The **AMT** is a high-power Class-IV laser system as described by the **Center for Devices and Radiological Health (CDRH)**.

Warning

The accessible laser radiation is very dangerous and can cause serious damage to your eyes and skin. Even diffuse scattered radiation can be dangerous.

Wear Laser Safety Glasses!

The laser radiation generated by this laser system is in the invisible range (1064 nm). If the system is equipped with the second harmonic module, third harmonic or forth harmonic then the system will also generate laser radiation in the visible range (532 nm) and ultra-violet (355 nm).

Nominal Ocular Hazard Area

Warning

Even apparently defused surfaces may reflect laser radiation.

Defuse reflected laser beams and multiple reflected laser beams are hazardous.

Always wear laser safety glasses when working on or with the **AMT**. Your safety glasses should be optically dense at 1064, 532, and 355nm. Please consult your Laser Safety Officer for the correct glasses that conform to your company's requirements.

The laser area must be clearly marked with appropriate warning signs and/or warning lamp or secured with laser safe blinds and/or interlocks.

Make sure there are no flammable or explosive objects and/or materials inside the nominal ocular hazard area.

Laser radiation may cause fires and explosions.

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Laser safety goggles

Always wear safety goggles when working on the StarCut Tube SL Femto, as wave lengths 900-1250 nm are emitted, [see 9.2 Safety equipment on page 9-1!](#)

Protective gloves

Wear protective gloves to prevent injuries caused by the unplanned crate boards, supplemental items and sharp tools!

Safety shoes

Wear safety shoes for protection against sharp objects, such as nails, sheet metal, etc.!


Hard hat


Wear a hard hat when loading and unloading the laser system to prevent head injuries caused by falling objects!

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
- ✓ Laser (Nano, Pico, Femto)
- ✓ Beam Delivery Part
 - Direct Mode
 - Scanner Mode
- ✓ Vision part
- ✓ Stage Part
- ✓ Power Control Part

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


- A Motion Stage
- B Scanner
- C Laser Source

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System Control Commands

COMMAND FROM PC (Syntax)	RESPONSE FROM CONTROLLER	COMMENTS
MCH	1	000001: remote mode
DI1=0	0	Diode off
DI1=1	1	Diode on
QSW=0	0	Q-switch off
QSW=1	1	Q-switch on
SHT=0	0	Shutter off
SHT=1	1	Shutter on
EXT=0	0	Internal trigger
EXT=1	1	External trigger
PRF=xxxxxx	xxxxxx	Pulse rate set up to 300KHz. Unit is Hz.
PRF = 12345	12345	Sets the repetition rate to 12345 Hz (shown as 12. 3kHz
PRF = 123456	123456	Sets the repetition rate to 123456 Hz (shown as 123 kHz
PWF=xxxxxx	xxxxxx	Set power control factor. Range 1-980 corresponding to Minimum-Maximum power output. (x = 0 to 9)
PWF=500	500	Sets laser power output to approximately 50% of Maximum power
PWF=980	950	Sets laser power output to approximately 100% of Maximum.
IS1=xxxxxx	xxxxxx	Diode current set. Unit: 0.1A. . (x = 0 - 9)
IS1=50	50	Set diode current to 5.0A
IS1=21	21	Set diode current to 2.1A
PWS=xxx	xxxxxx	Functions only with optional power monitor Feature. Power calibration setting. Unit: 1mW/0.1W/0.01W (Unit depends on Power Monitor Type.) (x = 0 to 9)
PWS=888	888	Set calibration power to 8.88W (Power monitor for laser power <10W)
T1S=xxxxxx	xxxxxx	Sets the temperature for Crystal 1
T1S=4751	4751	Sets the crystal 1 set temperature. Example 4751 = 47.51C
TS2=xxxxxx	xxxxxx	Sets the temperature for Crystal 2
TS2=5200	5200	Sets the crystal 2 set temperature. Example 5200 = 52.00C
GTR=0	0	Set the Gate logic to TTL low is lasing
GTR=1	1	Set the Gate logic to TTL high is lasing
MOD=xxxxxx	xxxxxx	Sets the mode of the laser
MOD=12	12	Set the laser operation mode 12

System Query Commands

QUERY FROM PC	RESPONSE FROM CONTROLLER	Comments
?DI1	0 1	Diode off Diode on
?QSW	0 1	Q-switch off Q-switch on
?SHT	0 1	Shutter Closed Shutter Open
?EXT	0 1	Internal trigger mode External trigger mode
?PRF Example:	xxxxxx 23456	Pulse rate set. Unit HZ. Current repetition rate is 23,456Hz or 23.4kHz
?IS1 Example:	xxxxxx 50	Diode current setting. Unit: 0.1A. Reports the ste diode current setting is 5.0 amps
?IDI Example:	xxxxxx 50	Actual diode current set. Unit: 0.1A. Reports the actual diode current is 5.0 amps
?PWF Example:	xxxxxx 950	Power control factor (Range: 1—980. 980 for highest laser power output) Reports that power control factor is set to 950 - please refered to your factory supplied test data to understand the power factor curve.
?AOW Example:	xxxxxx 123	Hour meter reading. Unit: 1 hour. Report the current operation hours
?OKK Example:	2731	Firmware version number Displays the current firmware version - 2.7.3.1
?TK1	xx4751	Query Crystal 1 Actual temperature in units of 0.01 degC. Example 4751 = 47.51C
?T1S	xx4751	Query Crystal 1 Set temperature in units of 0.01 degC. Example 4751 = 47.51C
?TK2	xx5000	Query Crystal 2 Actual temperature in units of 0.01 degC. Example 5000 = 50.00C
?T2S	xx5000	Query Crystal 2 Set temperature in units of 0.01 degC. Example 5000 = 50.00C
?TLA		Query Actual laser head temperature in units of 0.01 degC.

System Query Commands, Continued



QUERY FROM PC	RESPONSE FROM CONTROLLER	Comments
?TDI		Query Diode 1 Actual temperature in units of 0.01 degC.
?GTR		000000: gate low trigger 000001: gate high trigger
?MOD		Query the laser operation mode. Example 000012 = Mode 12
?STA	0	Normal
	1	Diode 1 Temperature High
	2	Diode 1 Temperature Low
	3	Laser Head Temperature High
	4	Laser Head Temperature Low
	5	Environmental Temperature High
	6	Environmental Temperature Low
	7	Interlock Open
	8	Fiber error
	9	Diode 2 Temperature High
	10	Diode 2 Temperature Low
	11	Diode 3 Temperature High
	12	Diode 3 Temperature Low
	13	Crystal 1 Temperature High
	14	Crystal 1 Temperature Low
	15	Crystal 2 Temperature High
	16	Crystal 2 Temperature Low
	17	Crystal 3 Temperature High
	18	Crystal 3 Temperature Low
	19	Crystal 4 Temperature High
	20	Crystal 4 Temperature Low
	21	Not used
	22	Not used
	23	Not used
	24	Emergency button
	25	Not used

Special Controller Messages

SPECIAL CASE	RESPONSE FROM CONTROLLER	COMMENTS
Any error message: wrong message, wrong message ending (0xD, 0xA), or incorrect value range.	ERROR	Need to check the message, the ending and value. If it is caused by RS232 bit error rate, PC needs to resend message.
PC tries to send a message to laser controller when it is busy with other high priority task.	BUSY	Wait for laser driver to finish its task.
When laser alarms any command sent to laser driver will be ignored and the system will send back "ALARM", followed by the alarm condition code.	Example: ALARM '000007	Indicates laser is alarming because of interlock open need to verify
When alarm condition is cleared, the laser will require user to enter the command RTS=1	RESET	RESET by user required

4) AOC PICO GUI(Graphical User Interface) Software

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
5) AOC PICO GUI(Graphical User Interface)

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
Chiller – Turn on



Step 1: Turn on the switch

Step 2: Press the power button
- Wait for about 30 minutes for the actual temperature to match the set temperature. (22.0 °C)

(2) Laser turn on- Laser head

Only turn on the Laser head after Chiller temperature has stabilized



Step 1: Press the power button
- Wait for 10 second

Step 2: Switch the key

Laser turn on – AOC software

(3) Turn on Power Supply

Step 1: Switch the power button

Step 2: Turn on the switch


All motion are ready to use.

(4) Turn on PC

Turn on PC Power

(5) AOC PICO GUI(Graphical User Interface)

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

1) KOS Software


2) Machining

1. Connect to Laser
2. Select the Alignment setting → Check Auto Focus
3. Select the Item number.
4. Input the item thickness.
5. The Connect to the attenuator.
6. Input a value and then press the “Enter” key and Shutter on
7. Open/Generate a CAD Drawing

3) Start machining.


1. Select Jog Motion or Relative Motion
2. IF: Relative Motion Distance of Camera X/Y/Z-axis and Stage Z-axis
3. Select Right/Left Camera for motion control.
4. Execute Camera X/Y/Z-axis Motion by clicking button.
5. Select Right/Left Stage for motion control.
6. Execute Stage Z-axis Motion by clicking button


4) Introduction of Machining Group

Machining Group

- Before machining, users can select one machining group for each DXF Object (etc Line, Circle, LwPolyline..).
- User can set the specified parameters in each machining group.
- Main Parameters :
 1. Processing Mode -> Optional
 2. Processing Speed (mm/s)
 3. Acceleration (mm/s²) -> Optional
 4. Attenuator Power Level (%) -> Optional
 5. Repetition Rate (kHz)
 6. Number of Pulses -> Optional
 7. MCS (LaserPower level)
 8. Control Signal -> Optional
 9. Spot Dia. (um) -> User need input it in advance
 10. LSO (%) (Laser Spot Overlap)
 11. Fluence (J/cm²)
 12. Power Energy (uj)
 13. Average Power (Watt)
 14. Machining Z Offset (mm)
 15. Pass #
 16. Pass X Offset (mm) -> Optional
 17. Pass Y Offset (mm) -> Optional
 18. Laser On Delay (ms) -> Optional
 19. Laser Off Delay (ms) -> Optional
 20. Group Delay (s)
 21. Processing Mode -> Optional

Machining Group – Common Operation

Operation Example

1. Input specified parameters in Group 3
2. Specify Line 5 to Line 10 for machining by using Group 3's machining parameters.
3. Update group number
4. Update group color
5. Select Machining Mode (Optional)
 1. On the Fly (IFOV Laser Machining)
 2. Only Scanner
 3. Only Stage
6. Perform Machining

Machining Group – Processing Speed (mm/s)

Alignment Tool (Scanner Mode)

- Processing Speed mean a scanner's scan speed during machining.

Processing Speed (mm/s)
150
150
150
150

Machining Group – Laser Parameters

- Repetition Rate is Laser Pulse Repetition Rate (PRF)
- Min./Max. Value depend on the specification of Laser.

Repetition Rate (kHz)
300
500

- MCS is Laser Power level (Power Control)
- Power change is non-linear.

MCS
300
1000

- Attenuator Power Level can meet the demand of sensitive film processing
- Power change is linear.
- Power Level Range : 0% ~ 100%

Attenuator Power Level(%)
20
35

- Beam Expansion and Divergence are used in MBE (Motorized Beam Expander) function
- The 2 parameters are optional.




B. Expansion	Divergence
20001	20002
20000	20000

Machining Group – Laser Parameters

- Number of Pulses is used in Burst Mode
- High level gate mode is the normal Control Signal
- The 2 parameters are optional.
- Using Rising edge trigger mode (Control Signal) and Number of Pulses can be used in TGV (Through Glass Vias)

Number of Pulses	Control Signal
1	High level gate mode
1	High level gate mode

Machining Group – Laser Parameters


LSO Calculation

- LSO (%) is Laser Spot Overlap
- The LSO formula is as below
- Input 3 variables of v,f,d to get LSO

Laser Spot Overlap

$$\text{LSO (\%)} = \left(1 - \frac{v}{f \times d}\right) \times 100\%$$


LSO = Laser Spot Overlap
 v = processing speed (mm/s)
 f = repetition frequency (kHz)
 d = beam spot diameter (um)


LSO = 50%

LSO = 0%

LSO = -50%


no.	Processing Speed (mm/s)	PRF (kHz)	Power Level	Spot Dia.(um)	LSO(%)
1	9600	400	300	20	-20

Input

Input

Input

Output

Fluence Calculation

- Fluence is optical energy delivered per unit area
- The fluence formula is the below
- Use three variables of p,f,d to change Fluence

Fluence

Fluence = 5.427056 J/cm²

$$\text{Fluence (J/cm}^2\text{)} = p / f / (\pi * (d / 2)^2)$$

p = laser power (Watt)

13.32 Watt

f = repetition frequency (kHz)

500 kHz

d = beam spot diameter (um)

25 um

Calculate

PRF (kHz)	Power Level	Spot Dia.(um)	LSO(%)	Power (Watt)	Fluence(J/cm ²)
400	300	20	-20	20	7.95774715459477

Input

Input


Input

Output

Machining Group – Machining Z Offset (mm)

For Alignment Tool (Scanner Mode)

- Machining Z Offset is Z-axis move during machining by controlling the Z-axis position of Sample stage.
- Users can get the optimal spot size by using this function.


Machining Group – Pass Parameters

- Pass # is the repeat times of machining.

Pass #	
3	1
2	
1	

The Pass # value of Machining Group 2 is 2.

Each line path of Machining Group 2 will be passed for 2 times by laser beam.

- Pass X Offset is the added Offset in every repeat times of machining in X direction.
- Pass Y Offset is the added Offset in every repeat times of machining in Y direction.
- The 2 parameters are optional.


Pass #	Pass X Offset (mm)	Pass Y Offset (mm)
3	0.2	0.6
2	0.1	0.5
1	0	0

Machining Group – Group Parameters

- Group Delay (s) is the delay time of start machining of each machining group.

Group Delay (s)
5
10
30

The group delay value of Machining Group 3 is 30s.


Wait for 30s before processing the line path of Machining Group 3.

- Group Note : User can record the purpose for each Machining Group.

no.	Note	Be
1	Wafer use	Gr
2	Glass use	IR
3	Soft substrate use	IR

Machining Group – Processing Mode


- Processing Mode have 3 kind of machining mode for DXF 2D Polyline object.
 - Outline
 - Filling
 - Outline and Filling


Filling Setting

Hatching


Pitch (mm) : 0.5 ☒ Cross Hatch
 Angle (°) : -10
 Offset (mm) : 0
 Apply

Outline

Outline and Filling

Filling

Machining Group – Export All Group Parameters

no.	Note	Beam Path	Processing Mode	Processing Speed (mm/s)	Attenuator Power Level(%)	Repetition Rate (kHz)	MCS	Number of Pulses	Control Signal
1	Wafer use	Green IFOV Convex	Outline and Filling	100	100	999	300	1	High level g...
2	Glass use	IR	Outline	100	0	500	1000	1	High level g...
3	Soft substrate use	IR	Outline	200	0	500	1000	1	High level g...
4		IR	Outline	100	0	500	300	1	High level g...
5		IR	Filling	100	0	500	300	1	High level g...

☒ Control Laser Parameter ☒ Control Attenuator Power ☒ Control MBE ☐ Control Acceleration Use Max. Acc : 10000mm/s^2

Export all parameter to a CSV file


A	B	C	D	E	F	G	H	I	J
no.	Note	Beam Path	Processing Mode	Processing Speed (mm/s)	Attenuator Power Level(%)	Repetition Rate (kHz)	MCS	Number of Pulses	Control Signal
1	1	Wafer use	100	100	999	300	20001	20002	20
2	2	Glass use	100	0	500	1000	20000	20000	20
3	3	Soft substrate use	200	0	500	1000	20000	20000	20
4	4		100	0	500	300	20000	20000	20
5	5		100	0	500	300	20000	20000	20
6	6		300	0	500	300	20000	20000	20
7	7		100	0	500	300	20000	20000	20

Open CSV file by Excel

5) Device Control

KOS software can intergrate with many device for customer's demand


Device Control

KOS Software can integrate with many devices for customer's demand.

1. Laser Control
2. Power Meter Control
3. Attenuator Control
4. Camera
5. Motorized Beam Expander (MBE) -> **Optional**
6. Laser Displacement Sensors (LDS) -> **Optional**
7. Machining Group can automatically control the above some devices function during machining.

Device Control - Laser Control

1. Select Com Port
2. Connect to Laser
3. Turn On Laser and then wait for 8 mins
4. Select Wavelength and then wait for 2 mins
5. Turn On / Off laser shutter
6. Turn On Sleep Mode if user doesn't use laser for a long time. It can extent laser use life.
7. Check the check box of Interpolat and then user can input the value of Repetition Rate.
8. Select the value of Repetition Rate
9. Input the value of Power Control (Change Power)
10. Input Number of Pulses/Burst -> **Optional**
11. Control Signal Setting
 - High level gate mode -> **Normal Use**
 - Low level gate mode
 - Rising edge trigger mode -> **Use it with Number of Pulse/Burst**
 - Falling edge trigger mode
12. Our software integrated with the above some laser control function
13. Laser Fire -> Beam Appear, Laser Stop -> Beam Disappear

Device Control – Common Device Operation

Item Setting Laser / Power Meter/ Attenuator / MBE / LDS Control

AOC Laser LDS Control

1 Connect

Shutter: 2 on off

Repetition Rate (kHz): 550

Number of Pulses / Burst: 3

MCS / Power Control: 100

Serial Port: None

Control Signal Setting:

Machining Group Filling Setting Laser

Power Meter Control

4 Connect

Average Power (Watt): 7

14 Disconnect

Ready Pos

5 Measure Pos

Wavelength (nm)

6 Laser Fire

Laser Stop 13

Attenuator Control

8 Connect

Power Level (%): 9 25

Wavelength (nm)

Custom Table:

Note	B. Expansion	Div
1.2x	-7787687	564
3.6x	-2343467	454
5x	-5656576	756

MBE Control

10 Connect

Beam Diameter: 1.5

Beam Expansion: 10

Divergence: 1.5


Common Device Operation

1. Connect to Laser
2. Turn on Laser Shutter
3. Check / Input Laser parameters
4. Connect to Power Meter
5. Move Power Meter to Measure Position
6. Laser Fire
7. Check Power Value
8. Connect to Attenuator
9. Input a value (%)
 - Enter key
 - Check Power Value
10. Connect to MBE
11. Input MBE parameters
 - Check Beam Shape
 - Check Power Value
12. Record MBE test result in MBE Table Tab
13. Laser Stop
14. Move Power Meter back to Ready Position


1) Scientific Chart

Scientific Chart

- KOS Software provide 3 scientific charts to meet customer demand.

Scientific Chart – LSO 3D Mesh Chart

LSO 3D Mesh Chart

• User can use it to get the optimal LSO conditions after program calculation

• Common Operation

1. Input Parameters
2. Generate Data List and 3D Chart
3. Edit Data
4. Update 3D Chart
5. Export Chart Table
6. Show Mesh
7. Export Mesh Table Data

Scientific Chart – Power, Fluence and Pulse Energy vs. Machining Group

• Users can check the actual condition of Power, Fluence and Pulse Energy for each laser parameters in Machining Group before machining.

• Common Operation

1. Connect to Laser
2. Turn on shutter
3. Connect to Power Meter
4. Move Power Meter to Measure Position
5. Edit laser parameters in Machining Group
6. Input Machining Group No. Range
7. Start Measure Power
8. Export Graph
9. Move Power Meter to Ready Position

Scientific Chart – Power and Pulse Energy vs. Repetition Rate

6) Measurement

Small Size Measure

- Small Size Measure : Use Camera to measure Distance, Angle, Area, ...

Direct Writing Tool

Alignment Point 1 and Alignment Point 2

Select AP1 Mark	Save AP1 to Get AP2	AP1	AP2	Select AP2 Mark	Click
X1 Offset:	AP1 X (mm): 316.902	AP2 X (mm): 416.902	X2 Offset:		
Y1 Offset:	AP1 Y (mm): -151.06	AP2 Y (mm): -151.06	Y2 Offset:	Measure	


Alignment Tool

Alignment Point 2

Select Mark	AP2 X (mm): -44.11	Center2 X (mm): -49.11
X2 Offset:	AP2 Y (mm): -69.297	Center2 Y (mm): -74.297
Y2 Offset:	Go to AP2	Center Offset: 40 mm
Calculate Approx. AP2	Save AP2	Measure
		Center2
		Save Center2

Alignment Point 1

Select Mark	AP1 X (mm): 26.081	Center1 X (mm): 31.081
X1 Offset:	AP1 Y (mm): -69.982	Center1 Y (mm): -73.982
Y1 Offset:	Go to AP1	Center Offset: 40 mm
Calculate Approx. AP1	Save AP1	Measure
		Center1
		Save Center1

Big Size Measure

Direct Writing Tool

- Use Linear / Gantry Stage and Camera to measure size, angle and so on.

Large Area Measuring / Angle Correction			
1 Save 1st Point	1st X (mm) : 35.874	1st Y (mm) : -31.184	1st Point
2 Save 2nd Point	2nd X (mm) : 243.439	2nd Y (mm) : -174.271	2nd Point
Start Angle Correction	Mid X (mm) : 235.872	Mid Y (mm) : -174.264	Mid Point
Recover Correct Angle	Angle (°) : -0.057	Dist. (mm) : 15.135	Auto-Correct
Save After Point	After X (mm) : 0	After Y (mm) : 0	Point
Get Rotation Center	Center X (mm) : 0	Center Y (mm) : 0	Point

Measure Operation

1. Move Camera to 1st Point and then save the position
2. Move Camera to 2nd Point and then save the position
3. Show the calculation result
4. Automatically Move Camera to the position of 1st Point after Angle Correction


Alignment Tool

- Use Linear / Gantry Stage and Camera to measure size, angle and so on.

Large Area Measuring & Angle Correction			
1 Save 1st Point	Angle (°)	Correct Angle	1st X (mm)
2 Save 2nd Point	Distance (mm)		1st Y (mm)
1st Point	Mid X (mm)	Mid Point	2nd X (mm)
2nd Point	Mid Y (mm)		2nd Y (mm)

Auto-Correct (Optional)

- Make a line horizontal with camera by stage rotation and then move camera to the 1st Point from a former 2nd Point.

7) CAD Generator

Machine Vision | Machining | Scientific Chart | Image to DXF | CMD

1. Click the icon of CAD Generator

2. Input the parameters according to user needs

3. Click Generate Button

4. Users can directly perform machining by using the drawing.

CAD Generator - Please input parameters

Basic Line Array			
Horizontal line number :	12		
Horizontal line length (mm) :	100		
Horizontal line interval (mm) :	10		
Horizontal group number :	0		
Horizontal group interval (mm) :	0		
Vertical line number :	0		
Vertical line length (mm) :	1		
Vertical line interval (mm) :	3.3		
Vertical group number :	0		
Vertical group interval (mm) :	0		
<input type="checkbox"/> Group interval include line interval			
<input type="checkbox"/> Hide Line Array			

Circle Array			
Horizontal Circle number :	1		
Horizontal Circle Radius (mm) :	1		
Vertical Circle number :	0		
Vertical Circle Radius (mm) :	1		
<input type="checkbox"/> Circle Size depend on Beam Size (Percussion Drilling)			
<input checked="" type="checkbox"/> Hide Circle Array			

Rectangle Array			
Horizontal Rectangle number :	1		
Horizontal Rectangle Width (mm) :	1		
Horizontal Rectangle Height (mm) :	1		
Horizontal Rectangle Corner R (mm) :	0.2		
Vertical Rectangle number :	0		
Vertical Rectangle Width (mm) :	10		
Vertical Rectangle Height (mm) :	10		
Vertical Rectangle Corner R (mm) :	0		
<input checked="" type="checkbox"/> Hide Rectangle Array			

Helical Circle Array			
Horizontal Helical Circle number :	1		
Horizontal Helical Circle Radius (mm) :	20		
Horizontal Helical Circle Interval (mm) :	1		
Vertical Helical Circle number :	0		
Vertical Helical Circle Radius (mm) :	10		
Vertical Helical Circle Interval (mm) :	1		
<input checked="" type="checkbox"/> Hide Helical Circle Array			

Total Setting

Default Machining group No : 1

Generate Cancel

Item Setting | Laser / Power Meter / Attenuator / MBE / LDS Control | Machining

no.	Note	Beam Path	Processing Mode	Processing Speed
1	Wafer use	Green IFOV Convex	Outline and Filling	100
2	Glass use	IR	Outline	100
3	Soft substrate use	IR	Outline	200
4		IR	Outline	100
5		IR	Filling	100

Crosshair Size (mm) : 0.5 Position Drawing of AP1 & AP2

CAD Generator

Automatically generate Line Array, Circle Array, Rectangle Array and Helical Circle Array by using the parameters that user input.

Basic Line Array, Circle Array, Rectangle Array and Helical Circle Array can be generated at the same

Common Operation:

1. Click the icon of CAD Generator
2. Input the parameters according to user needs
3. Click Generate Button
4. Users can directly perform machining by using the drawing.

CAD Generator – Basic Line Array


Basic Line Array	
Horizontal line number :	12
Horizontal line length (mm) :	100
Horizontal line interval (mm) :	10
Horizontal group line number :	3
Horizontal group interval (mm) :	50
Vertical line number :	0
Vertical line length (mm) :	0
Vertical line interval (mm) :	0
Vertical group line number :	0
Vertical group interval (mm) :	0
<input type="checkbox"/> Group interval include line interval	
<input type="checkbox"/> Hide Line Array	

Basic Line Array	
Horizontal line number :	12
Horizontal line length (mm) :	100
Horizontal line interval (mm) :	10
Horizontal group line number :	3
Horizontal group interval (mm) :	50
Vertical line number :	0
Vertical line length (mm) :	0
Vertical line interval (mm) :	0
Vertical group line number :	0
Vertical group interval (mm) :	0
<input checked="" type="checkbox"/> Group interval include line interval	
<input type="checkbox"/> Hide Line Array	

Basic Line Array	
Horizontal line number :	0
Horizontal line length (mm) :	0
Horizontal line interval (mm) :	0
Horizontal group line number :	0
Horizontal group interval (mm) :	0
Vertical line number :	12
Vertical line length (mm) :	100
Vertical line interval (mm) :	10
Vertical group line number :	3
Vertical group interval (mm) :	50
<input checked="" type="checkbox"/> Group interval include line interval	
<input type="checkbox"/> Hide Line Array	

Basic Line Array

- Horizontal Lines
- Vertical Lines
- Circle Array, Rectangle Array and Helical Circle Array follow Basic Line Array Setting to generate Circles, Rectangles and Helical Circles and then hide Basic Line Array.


CAD Generator – Circle Array

Basic Line Array	
Horizontal line number :	12
Horizontal line length (mm) :	100
Horizontal line interval (mm) :	10
Horizontal group line number :	3
Horizontal group interval (mm) :	50
Vertical line number :	0
Vertical line length (mm) :	0
Vertical line interval (mm) :	0
Vertical group line number :	0
Vertical group interval (mm) :	0
<input type="checkbox"/> Group interval include line interval	
<input checked="" type="checkbox"/> Hide Line Array	

Circle Array	
Horizontal Circle number :	3
Horizontal Circle Radius (mm) :	4
Vertical Circle number :	0
Vertical Circle Radius (mm) :	0
<input type="checkbox"/> Circle Size depend on Beam Size (Percussion Drilling)	
<input type="checkbox"/> Hide Circle Array	

Total Setting	
Default Machining group No. :	1
Generate Cancel	

Circle Array

- Horizontal Circles
- Vertical Circles
- Common Operation
 1. Input Basic Line Array Setting
 2. Hide Line Array
 3. Input Circle Array Setting
 4. Show Circle Array
 5. Input Machining group No.
 6. Generate Circle Array

- Circle Size depends on Beam Size
 - TGV application
 - This function is optional.
 - This function must be use with the setting of Number of Pulses and Control Signal

Number of Pulses	Control Signal
2	Rising edge trigger mode


Pulse Duration : 207 fs

1 shot 2 shot 3 shot 4 shot

Pulse Count Count Count Count


= 1 = 2 = 3 = 4

CAD Generator – Rectangle Array


Rectangle Array

- Horizontal Rectangles
- Vertical Rectangles
- Common Operation
 1. Input Basic Line Array Setting
 2. Hide Line Array
 3. Input Rectangle Array Setting
 4. Show Rectangle Array
 5. Input Machining group No.
 6. Generate Rectangle Array


CAD Generator – Rectangle Array

Rectangle Array

- Horizontal Rectangles
- Vertical Rectangles
- Common Operation
 1. Input Basic Line Array Setting
 2. Hide Line Array
 3. Input Rectangle Array Setting
 4. Show Rectangle Array
 5. Input Machining group No.
 6. Generate Rectangle Array

CAD Generator – Helical Circle Array

Basic Line Array

Horizontal line number :	6
Horizontal line length (mm) :	100
Horizontal line interval (mm) :	30
Horizontal group line number :	3
Horizontal group interval (mm) :	50
Vertical line number :	0
Vertical line length (mm) :	0
Vertical line interval (mm) :	0
Vertical group line number :	0
Vertical group interval (mm) :	0

☐ Group interval include line interval
☒ Hide Line Array

Helical Circle Array

Horizontal Helical Circle number :	6
Horizontal Helical Circle Radius (mm) :	4
Horizontal Helical Circle Interval (mm) :	2
Vertical Helical Circle number :	0
Vertical Helical Circle Radius (mm) :	10
Vertical Helical Circle Interval (mm) :	1

☐ Hide Helical Circle Array

Total Setting

Default Machining group No. : 4

Generate Cancel

Horizontal Helical Circle Number : 6
 -> It means generate 4 helical circles on each horizontal line.

Horizontal Helical Circle Number : 6
 -> It means generate 4 helical circles on each horizontal line.


- ## Helical Circle Array
- Horizontal Helical Circles
 - Vertical Helical Circles
 - Common Operation
 1. Input Basic Line Array Setting
 2. Hide Line Array
 3. Input Helical Circle Array Setting
 4. Show Helical Circle Array
 5. Input Machining group No.
 6. Generate Helical Circle Array

7) Image to DXF


Image to DXF

- KOS Software provide the function to convert an image file to DXF file.
- Common Operation :
1. Select an image file (bmp, png, jpg).
 2. Edit Conversion Parameters.
 3. Convert Image to Binarization Image. If the image result is not good , try Step. 2 again.
 4. Convert the image to DXF.
 5. Perform an image machining .


8) Alignment Tool – Motion Control

1. Jog Motion is when motion is commanded at a defined constant velocity with infinite distance.
2. Select to control a Left or Right Camera Stage
3. Move Camera Stage by pressing buttons
4. Select UV/Green/IR Sample Stage if there is not only one sample stage
5. Move Sample Stage by pressing buttons

1. Relative motion is a motion to move by a specified distance from the current position.
2. Input a specified relative distance of X/Y/Z/StageZ axis
 - o X/Y/Z : Camera Stage Axis
 - o Stage Z : Sample Stage Axis
3. Select to control a Left or Right Camera Stage
4. Move Camera Stage by clicking buttons
5. Select UV/Green/IR Sample Stage to control motion if there is not only one sample stage
6. Move Sample Stage by clicking buttons


Alignment Tool :


1. Alignment by using 2 Camera
2. Use Left Camera to move to the position of Alignment Point 1 (AP1)
3. Use Right Camera to move to the position of Alignment Point 2 (AP2)
4. It's for the Beam Path of Scanner Mode
5. Mark Alignment
6. Edge Alignment

Alignment Tool

Definition – Aligning Width, Alignment Height


Alignment Tool – Initial Setting


1. Execute Homing
2. Select Item
3. Select Beam Path
4. Input Aligning Width
5. Input Aligning Height
6. Input Item Thickness

Alignment Tool – P1 Processing

1. New Sample Init
2. Uncheck / Check Edge Alignment
3. Move Left Camera to AP1
4. Move Right Camera to AP2
5. Put sample marks/edges under Camera Live views

Alignment Tool – P1 Processing

Alignment Tool – P1 Processing

Alignment Point 2

Select Mark

1. Set : 0 / 0

Y2 Offset : 0 / 0

Calculate Approx. AP2

Mark Alignment

Edge Alignment

1. Select Mark
2. Select a mark center/edge by using mouse or "up", "down", "left" and "right" keys
3. Get Offsets

Alignment Tool – P1 Processing

File Path : C:\KOS\code\AlignmentTool\bin\Debug\GenLines\genLines-22-0-28.dxf



Object no.	Type	Machining Group
1	Line	3
2	Line	3
3	Line	3

Aligning Item	Aligning Height (mm)	Item Thickness (mm)	AP1 X (mm)	AP1 Y (mm)	AP2 X (mm)	AP2 Y (mm)	User X Offset (mm)	User Y Offset (mm)	User Angle Offset (°)	Auto Mach
1	100	1.1	25.561	-58.962	-44.11	49.3297	0	0	22.653	
2	150	0.8	0	0	0	0	0	0	-2.653	
3	150	1.3	0	0	0	0	0	0	-0.5	
4	0	0	0	0	0	0	0	0	0	

1. Open/Generate a CAD
2. Input User X Offset
3. Input User Y Offset
4. Start machining. After machining, the user can take out the sample.

Alignment Tool – P1 Processing

P1 Processing Result

Alignment Tool – P2/P3 Processing

1. Put Marks /Edges under Camera Live Views

Alignment Tool – P2/P3 Processing

Alignment Point 1

Select Mark
 1 set :
 Y1 Offset :
 Calculate Approx. AP1

1. Select Mark
 2. Program select the former position automatically
 3. Get Offsets

Mark Alignment

Edge Alignment

Alignment Tool – P2/P3 Processing

Alignment Point 2


Select Mark
 1 fset :
 Y2 Offset :
 Calculate Approx. AP2

1. Select Mark
 2. Program select the former position automatically
 3. Get Offsets

Mark Alignment

Edge Alignment



Alignment Tool – P2/P3 Processing

1. Input User X Offset
2. Input User Y Offset
3. Start machining. After machining, the user can take out the sample.

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Error Code	Error Message	Interpretation
000001	!!! WARNING !!! TOO HIGH PUMP DIODE TEMP	Diode Temperature is higher than 35 °C
000002	!!! WARNING !!! TOO LOW PUMP DIODE TEMP	Diode Temperature is lower than 15 °C; or chiller controlling cable is not connected; or chiller water flow is too slow
000003	!!! WARNING !!! TOO HIGH LASER HEAD TEMP	Laser Head Temperature is higher than 40 °C
000004	!!! WARNING !!! TOO LOW LASER HEAD TEMP	Laser Head Temperature is lower than 15 °C
000005	!!! WARNING !!! TOO HIGH ENVIRONMENTAL TEMP	Environmental Temperature is higher than the setting alarm temperature (usually set to 32 °C or 35 °C)
000006	!!! WARNING !!! TOO LOW ENVIRONMENTAL TEMP	Environmental Temperature is lower than 15 °C
000007	!!! WARNING !!! INTERLOCK OPEN	Interlock is open

(3) LASER Controller Check

COMMAND FROM PC	RESPONSE FROM LASER DRIVER IS ALWAYS 6 DIGITS (EACH DIGIT IS BETWEEN 0-9)	COMMENTS
MCH	MOD_CH REMOTE LOCAL ERROR	If laser driver is in Local mode and LCD display stays at status screen, this command will change laser to REMOTE mode. Indicates laser driver is in Local mode but not at status screen Indicates laser driver is in REMOTE mode
MCH=1	MOD_CH LOCAL LOCAL	If laser driver is in REMOTE mode, this command will change laser driver to Local mode Indicates laser driver is in Local mode
DIO=0	000000	Diode off
DIO=1	000001	Diode on
QSW=0	000000	Q-switch off
QSW=1	000001	Q-switch on
SHT=0	000000	Shutter off
SHT=1	000001	Shutter on
EXT=0	000000	Internal trigger
EXT=1	000001	External trigger
PRF=xxxxxx (x: digits 0-9) Examples: PRF=12345 PRF=123456	xxxxxx 012345 123456	Pulse rate set up to 300KHz. Unit is HZ. Set pulse rate to 12345HZ (shown on display as 12.3KHZ) Set pulse rate to 123456HZ (shown on display as 123KHZ)
PWS=xxx (x: digits 0-9) Example: PWS=888	xxxxxx 000888	Power calibration setting. Unit: 1mW/0.1W/0.01W (Unit depends on Power Monitor Type.) Set calibration power to 8.88W (Power monitor for laser power <10W)
IDI=xxx (x: digits 0-9) Examples: IDI=150 IDI=321	xxxxxx 000150 000321	Diode current set. Unit: 0.1A. Set diode current to 15.0A Set diode current to 32.1A
PWF=xxxx (x: digits 0-9.) (Set power control factor. Range 1-1000 corresponding to Minimum-Maximum power output.) Examples: PWF=1000 PWF=950	xxxxxx 001000 000950	 Set laser power output to Maximum Set laser power output to about 95% Maximum.

QUERY FROM PC	RESPONSE FROM LASER DRIVER IS ALWAYS 6 DIGITS (EACH DIGIT IS BETWEEN 0-9)	COMMENTS
?DIO	000000 000001	Diode off Diode on
?QSW	000000 000001	Q-switch off Q-switch on
?SHT	000000 000001	Shutter off Shutter on
?EXT	000000 000001	Use internal pulse Use external pulse
?PRF	xxxxxx	Pulse rate set. Unit HZ.
Example: ?PRF	023456	Pulse rate is 23456HZ (23.4 KHZ)
?IDS	xxxxxx	Diode current setting. Unit: 0.1A.
Example: ?IDS	000150	Diode current setting is 15.0A
?IDI	xxxxxx	Actual diode current set. Unit: 0.1A.
Example: ?IDI	000150	Actual diode current is 15.0A
?PWS	xxxxxx	Power calibration setting. Unit: 1mW/0.1W/0.01W (Unit depends on Power Monitor Type.)
Example: ?PWS	000095	Power calibration setting: 9.5W (power monitor for laser power <10W)
?PWR	xxxxxx	Calibrated power measurement. Unit: 1mW/0.1W/0.01W (Unit depends on Power Monitor Type.)
Example: ?PWR	000095	Power monitor measurement : 9.5W (power monitor for laser power > 10W)
?PWF	xxxxxx	Power control factor (Range: 1—1000. 1000 for highest laser power output)
Example: ?PWF	001000	Laser output set to Maximum
?AOW	xxxxxx	Hour meter reading. Unit: 1 hour.
Example: ?AOW	000123	Hour meter reading is 123 hour.
?STA	000000 000001 000002 000003 000004 000005 000006 000007	Normal Diode Temperature High Diode Temperature Low Laser Head Temperature High Laser Head Temperature Low Environmental Temperature High Environmental Temperature Low Interlock Open

(4) System Control Commands

COMMAND FROM PC (Syntax)	RESPONSE FROM CONTROLLER	COMMENTS
MCH	1	000001: remote mode
DI1=0	0	Diode off
DI1=1	1	Diode on
QSW=0	0	Q-switch off
QSW=1	1	Q-switch on
SHT=0	0	Shutter off
SHT=1	1	Shutter on
EXT=0	0	Internal trigger
EXT=1	1	External trigger
PRF=xxxxxx PRF = 12345 PRF = 123456	xxxxxx 12345 123456	Pulse rate set up to 300KHz. Unit is Hz. Sets the repetition rate to 12345 Hz (shown as 12. 3kHz Sets the repetition rate to 123456 Hz (shown as 123 kHz
PWF=xxxxxx PWF=500 PWF=980	xxxxxx 500 950	Set power control factor. Range 1-980 corresponding to Minimum-Maximum power output. (x = 0 to 9) Sets laser power output to approximately 50% of Maximum power Sets laser power output to approximately 100% of Maximum.
IS1=xxxxxx IS1=50 IS1=21	xxxxxx 50 21	Diode current set. Unit: 0.1A. (x = 0 - 9) Set diode current to 5.0A Set diode current to 2.1A
PWS=xxx PWS=888	xxxxxx 888	Functions only with optional power monitor Feature. Power calibration setting. Unit: 1mW/0.1W/0.01W (Unit depends on Power Monitor Type.) (x = 0 to 9) Set calibration power to 8.88W (Power monitor for laser power <10W)
T1S=xxxxxx T1S=4751	xxxxxx 4751	Sets the temperature for Crystal 1 Sets the crystal 1 set temperature. Example 4751 = 47.51C
TS2=xxxxxx TS2=5200	xxxxxx 5200	Sets the temperature for Crystal 2 Sets the crystal 2 set temperature. Example 5200 = 52.00C
GTR=0 GTR=1	0 1	Set the Gate logic to TTL low is lasing Set the Gate logic to TTL high is lasing
MOD=xxxxxx MOD=12	xxxxxx 12	Sets the mode of the laser Set the laser operation mode 12

QUERY FROM PC	RESPONSE FROM CONTROLLER	Comments
?DI1	0 1	Diode off Diode on
?QSW	0 1	Q-switch off Q-switch on
?SHT	0 1	Shutter Closed Shutter Open
?EXT	0 1	Internal trigger mode External trigger mode
?PRF Example:	xxxxxx 23456	Pulse rate set. Unit HZ. Current repetition rate is 23,456Hz or 23.4kHz
?IS1 Example:	xxxxxx 50	Diode current setting. Unit: 0.1A. Reports the ste diode current setting is 5.0 amps
?DI Example:	xxxxxx 50	Actual diode current set. Unit: 0.1A. Reports the actual diode current is 5.0 amps
?PWF Example:	xxxxxx 950	Power control factor (Range: 1—980. 980 for highest laser power output) Reports that power control factor is set to 950 - please refered to your factory supplied test data to understand the power factor curve.
?AOW Example:	xxxxxx 123	Hour meter reading. Unit: 1 hour. Report the current operation hours
?OKK Example:	2731	Firmware version number Displays the current firmware version - 2.7.3.1
?TK1	xx4751	Query Crystal 1 Actual temperature in units of 0.01 degC. Example 4751 = 47.51C
?T1S	xx4751	Query Crystal 1 Set temperature in units of 0.01 degC. Example 4751 = 47.51C
?TK2	xx5000	Query Crystal 2 Actual temperature in units of 0.01 degC. Example 5000 = 50.00C
?T2S	xx5000	Query Crystal 2 Set temperature in units of 0.01 degC. Example 5000 = 50.00C
?TLA		Query Actual laser head temperature in units of 0.01 degC.

(5) System Query Commands, Continued

QUERY FROM PC	RESPONSE FROM CONTROLLER	Comments
?TDI		Query Diode 1 Actual temperature in units of 0.01 degC.
?GTR		000000: gate low trigger 000001: gate high trigger
?MOD		Query the laser operation mode. Example 000012 = Mode 12
?STA	0	Normal
	1	Diode 1 Temperature High
	2	Diode 1 Temperature Low
	3	Laser Head Temperature High
	4	Laser Head Temperature Low
	5	Environmental Temperature High
	6	Environmental Temperature Low
	7	Interlock Open
	8	Fiber error
	9	Diode 2 Temperature High
	10	Diode 2 Temperature Low
	11	Diode 3 Temperature High
	12	Diode 3 Temperature Low
	13	Crystal 1 Temperature High
	14	Crystal 1 Temperature Low
	15	Crystal 2 Temperature High
	16	Crystal 2 Temperature Low
	17	Crystal 3 Temperature High
	18	Crystal 3 Temperature Low
	19	Crystal 4 Temperature High
	20	Crystal 4 Temperature Low
	21	Not used
	22	Not used
	23	Not used
	24	Emergency button
	25	Not used

(6) Special Controller Messages

SPECIAL CASE	RESPONSE FROM CONTROLLER	COMMENTS
Any error message: wrong message, wrong message ending (0xD, 0xA), or incorrect value range.	ERROR	Need to check the message, the ending and value. If it is caused by RS232 bit error rate, PC needs to resend message.
PC tries to send a message to laser controller when it is busy with other high priority task.	BUSY	Wait for laser driver to finish its task.
When laser alarms any command sent to laser driver will be ignored and the system will send back "ALARM", followed by the alarm condition code.	Example: ALARM `000007	Indicates laser is alarming because of interlock open need to verify
When alarm condition is cleared, the laser will require user to enter the command RTS=1	RESET	RESET by user required

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